

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

| | | | |
|-----------------------------|---------------------------------|--------------------------------|---------------------------|
| Алматы (7273)495-231 | Казань (843)206-01-48 | Новокузнецк (3843)20-46-81 | Смоленск (4812)29-41-54 |
| Архангельск (8182)63-90-72 | Калининград (4012)72-03-81 | Новосибирск (383)227-86-73 | Сочи (862)225-72-31 |
| Астрахань (8512)99-46-04 | Калуга (4842)92-23-67 | Омск (3812)21-46-40 | Ставрополь (8652)20-65-13 |
| Барнаул (3852)73-04-60 | Кемерово (3842)65-04-62 | Орел (4862)44-53-42 | Сургут (3462)77-98-35 |
| Белгород (4722)40-23-64 | Киров (8332)68-02-04 | Оренбург (3532)37-68-04 | Тверь (4822)63-31-35 |
| Брянск (4832)59-03-52 | Краснодар (861)203-40-90 | Пенза (8412)22-31-16 | Томск (3822)98-41-53 |
| Владивосток (423)249-28-31 | Красноярск (391)204-63-61 | Пермь (342)205-81-47 | Тула (4872)74-02-29 |
| Волгоград (844)278-03-48 | Курск (4712)77-13-04 | Ростов-на-Дону (863)308-18-15 | Тюмень (3452)66-21-18 |
| Вологда (8172)26-41-59 | Липецк (4742)52-20-81 | Рязань (4912)46-61-64 | Ульяновск (8422)24-23-59 |
| Воронеж (473)204-51-73 | Магнитогорск (3519)55-03-13 | Самара (846)206-03-16 | Уфа (347)229-48-12 |
| Екатеринбург (343)384-55-89 | Москва (495)268-04-70 | Санкт-Петербург (812)309-46-40 | Хабаровск (4212)92-98-04 |
| Иваново (4932)77-34-06 | Мурманск (8152)59-64-93 | Саратов (845)249-38-78 | Челябинск (351)202-03-61 |
| Ижевск (3412)26-03-58 | Набережные Челны (8552)20-53-41 | Севастополь (8692)22-31-93 | Череповец (8202)49-02-64 |
| Иркутск (395)279-98-46 | Нижний Новгород (831)429-08-12 | Симферополь (3652)67-13-56 | Ярославль (4852)69-52-93 |
| Россия (495)268-04-70 | Киргизия (996)312-96-26-47 | Казахстан (7172)727-132 | |

Sonden – Probes
 Spatel – Watteträger
 Spatulas – Cotton Applicators

BN



1/2

BN184R

Spatel
 Spatula
 145 mm, 5 3/4"



16 mm

1/2

OM202R

Spatel
 Spatula
 145 mm, 5 3/4"

Ø 0,7 mm
 Diam. 0.7 mm



1/1

1/2

MODELL PICCOLO

OF261R

rund gerieft
 round serrated end
 110 mm, 4 1/2"

Ø 0,7 mm
 Diam. 0.7 mm



1/1

1/2

MODELL PICCOLO

OF262R

dreikant gerieft
 triangular serrated end
 110 mm, 4 1/2"

Ø 1,2 mm
 Diam. 1.2 mm



1/1

1/2

FARRELL

OF257R

rund gerieft
 round serrated end
 110 mm, 4 1/2"

Ø 1,4 mm
 Diam. 1.4 mm



1/1

1/2

LATHBURY-FARRELL

OF259R


dreikant gerieft
 triangular serrated end
 110 mm, 4 1/2"



Sonden – Probes
Watteträger
Cotton Applicators

BN

Ø 0,9 mm
Diam. 0.9 mm




BUCK

BN400R

spiralförmig gerieft
rough, spiral end
120 mm, 4 3/4"

Ø 0,9 mm
Diam. 0.9 mm



LATHBURY-FARRELL

BN467R


dreikant gerieft
triangular serrated end
140 mm, 5 1/2"



KINGS' COLLEGE HOSPITAL

BN403

spiralförmig gerieft
rough, spiral end
90 mm, 3 3/8"



KINGS' COLLEGE HOSPITAL

BN404

spiralförmig gerieft
rough, spiral end
135 mm, 5 1/8"

Ø 1,0 mm
Diam. 1.0 mm



BN440R

spiralförmig gerieft
rough, spiral end
140 mm, 5 1/2"

Ø 1,2 mm
Diam. 1.2 mm



BN445R

spiralförmig gerieft
rough, spiral end
140 mm, 5 1/2"




MILLIMETERS

Sonden – Probes
Watteträger
Cotton Applicators

BN



Ø 1,2 mm
Diam. 1.2 mm




FARRELL

BN460R

rund gerieft
round serrated end
140 mm, 5 1/2"

Ø 1,4 mm
Diam. 1.4 mm




LATHBURY-FARRELL

BN465R

dreikant gerieft
triangular serrated end
140 mm, 5 1/2"

Ø 1,2 mm
Diam. 1.2 mm



FARRELL

BN461R

spiralförmig gerieft
rough, spiral end
160 mm, 6 1/4"

Ø 1,2 mm
Diam. 1.2 mm




LATHBURY-FARRELL

BN466R

dreikant gerieft
triangular serrated end
160 mm, 6 1/4"

Ø 1,4 mm
Diam. 1.4 mm



LATHBURY-FARRELL

BN468R

dreikant gerieft
triangular serrated end
160 mm, 6 1/4"



Sonden – Probes
Watteträger
Cotton Applicators

BN

Ø 1,5 mm
Diam. 1.5 mm



BN452R

spiralförmig gerieft
rough, spiral end
180 mm, 7"


Ø 1,8 mm
Diam. 1.8 mm



BN457R

spiralförmig gerieft
rough, spiral end
180 mm, 7"


Ø 2,3 mm
Diam. 2.3 mm



HARTMANN

OM487R

spiralförmig gerieft
rough, spiral end
220 mm, 8 1/2"




MILLIMETERS



Wundhaken – Wound Retractors

Wundhaken, Lungenspatel
Retractors, Lung Spatulas

BT



SAUERBRUCH

BT514R

260 mm, 10 $1/4$ "

CORYLLOS

BT521R–BT523R

255 mm, 10"



MILLIMETERS

Wundhaken – Wound Retractors

Bauch- und Darmspatel

Abdominal and Intestinal Spatulas

BT



REVERDIN

BT749R

285 mm, 11 1/4"



HABERER

BT758R

biegsam
malleable
305 mm, 12"



RIBBON

| | |
|---------------|-------|
| BT756R | 25 mm |
| BT761R | 30 mm |
| BT762R | 40 mm |
| BT763R | 50 mm |

biegsam
malleable
330 mm, 13 1/8"



TUFFIER

BT750R

biegsam
malleable
200 mm, 8"



| | | | |
|---------------|-------|--------|-----|
| BT751R | 12 mm | 200 mm | 8" |
| BT752R | 17 mm | 200 mm | 8" |
| BT755R | 25 mm | 250 mm | 10" |

biegsam
malleable



Wundhaken – Wound Retractors
 Bauch- und Darmspatel
 Abdominal- and Intestinal Spatulas

BT




KADER

BT746R 30 mm

BT747R 40 mm

270 mm, 10 3/4"



biegsam
malleable

BT773R 30 mm

BT774R 40 mm

BT775R 50 mm

BT776R 60 mm

BT777R 70 mm

410 mm, 16"

MILLIMETERS



Wundhaken – Wound Retractors
Kropflöffel, Blutschöpfer
Depressor, Blood Scoops

BT




KOCHER

BT795R

Kropflöffel

Depressor, spoon shaped

220 mm, 8 3/4"



Plastic Surgery

Hooks and Retractors



BT966R

1/2



1/1



BT966R

MANNERFELT Double Spatula

| Catalog Number | Length (in.) |
|----------------|--------------|
| BT966R | 6 1/4 |



PLASTIC SURGERY




MILLIMETERS


Scheren – Scissors

DUROTIP® Ligaturscheren mit **W**ellenschliff

DUROTIP® **W**avecut Ligature Scissors



BC



speziell zum Schneiden von Nahtmaterial


- Hartmetalleinlagen gewähren dauerhafte Schnittschärfe (auch im distalen Bereich der Scherenblätter)
- der Spezial-**W**ellenschliff verhindert das Abgleiten des Nahtmaterials
- DUROTIP®-Ligaturscheren sind gekennzeichnet mit 1 vergoldeten Ring

especially for cutting suture materials

- Carbide inserts guarantee permanent cutting edge (also in the distal area of the blades)
- **W**avecut prevents suture materials from slipping out
- DUROTIP® ligature scissors are identified by 1 goldplated ring




BC294W-BC298W



BC294W

145 mm


5 3/4"



BC295W

180 mm


7"



BC296W

230 mm

9"



BC297W

260 mm

10 1/4"

BC298W

350 mm

14"

1/1

Klemmen – Forceps

Präparier- und Ligaturklemmen

Dissecting- and Ligature Forceps

BJ



BJ002R

Ligaturklemme gerade
Ligature forceps straight
95 mm, 3 3/4"



MICRO-ADSON

BJ008R

140 mm, 5 1/2"



BABY-OVERHOLT

BJ009R

135 mm, 5 1/4"



BABY-MIXTER

BJ011R-BJ012R

Obstetrics/Gynecology
Dilators / Depressors / Probes




BRAUN Uterine Depressor, Double Ended

| Catalog Number | Blade Size (mm) | Length (in.) |
|----------------|-----------------|--------------|
| E0030R | 24/28 | 10 5/8 |




E0030R

1/2



ENT-ORAL

Mouth Gags



OM171R-OM178R

RUSSEL-DAVIS Tongue Blade

| Catalog Number | Blade Size (mm) |
|----------------|-----------------|
| OM171R | 29 x 67 |
| OM172R | 33 x 75 |
| OM173R | 38 x 85 |
| OM174R | 40 x 92 |
| OM175R | 46 x 100 |
| OM176R | 18 x 57 |
| OM177R | 22 x 67 |
| OM178R | 25 x 75 |



OM171R



OM172R



OM173R



OM174R



OM175R



OM176R



OM177R







OM178R

ENT
ORAL

MILLIMETERS

General Surgery

Delicate Scissors



RAGNELL Dissecting Scissors

| Catalog Number | Tip | Blade | Length (in.) |
|----------------|-----|-------|--------------|
| BC026R | B/B | STR | 5 1/4 |
| BC027R | B/B | CVD | 5 1/4 |
| BC029R | B/B | CVD | 6 |



**OM222R**

LACK Tongue Depressor, curved, width: 22 mm, non-sterile, reusable

| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | LACK |
| Aesculap Product | |
| Groups | General Surgical Instruments, ENT |
| Geometry | curved |
| Width (mm) | 22.000 mm |
| Discipline | Maxillo-Facial Surgery, Ear-Nose-Throat Surgery |
| Indication | Tonsillectomy, Adenoidectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | Yes |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |

**OM640R**

MORITZ-SCHMIDT Tongue Depressor, 190 mm (7 1/2\"/>

| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | MORITZ-SCHMIDT |
| Aesculap Product | |
| Groups | ENT, General Surgical Instruments |
| Length | 190.00 mm |
| Length (Inch) | 7 1/2 |
| Discipline | Ear-Nose-Throat Surgery, Maxillo-Facial Surgery |
| Indication | Tonsillectomy, Adenoidectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | Yes |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |

**OM221R**


LACK Tongue Depressor, curved, width: 19 mm, non-sterile, reusable

| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | LACK |
| Aesculap Product | |
| Groups | General Surgical Instruments, ENT |
| Geometry | curved |
| Width (mm) | 19.000 mm |
| Discipline | Maxillo-Facial Surgery, Ear-Nose-Throat Surgery |
| Indication | Tonsillectomy, Adenoidectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | Yes |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |

**OM220R**


LACK Tongue Depressor, curved, width: 13 mm, non-sterile, reusable

| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | LACK |
| Aesculap Product | |
| Groups | General Surgical Instruments, ENT |
| Geometry | curved |
| Width (mm) | 13.000 mm |
| Discipline | Maxillo-Facial Surgery, Ear-Nose-Throat Surgery |
| Indication | Adenoidectomy, Tonsillectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | Yes |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |



OM212R

FRENZEL Tongue Depressor, 165 mm (6 1/2"), non-sterile, reusable



| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | FRENZEL |
| Aesculap Product | |
| Groups | General Surgical Instruments, ENT |
| Length | 165.00 mm |
| Length (Inch) | 6 1/2 |
| Discipline | Ear-Nose-Throat Surgery, Maxillo-Facial Surgery |
| Indication | Tonsillectomy, Adenoidectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | No |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |



OM236R

WEDER Tongue Depressor, 125 mm (5"), small, fenestrated, non-sterile, reusable



| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | WEDER |
| Aesculap Product | |
| Groups | ENT, General Surgical Instruments |
| Shape / Form | small |
| Length | 125.00 mm |
| Length (Inch) | 5 |
| Blades / Jaw | fenestrated |
| Discipline | Ear-Nose-Throat Surgery, Maxillo-Facial Surgery |
| Indication | Adenoidectomy, Tonsillectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | Yes |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |

**OM237R**

WEDER Tongue Depressor, 125 mm (5"), large, fenestrated, non-sterile, reusable



| | |
|-----------------------|---|
| Product search type | Tongue Depressor |
| Author | WEDER |
| Aesculap Product | |
| Groups | General Surgical Instruments, ENT |
| Shape / Form | large |
| Length | 125.00 mm |
| Length (Inch) | 5 |
| Blades / Jaw | fenestrated |
| Discipline | Maxillo-Facial Surgery, Ear-Nose-Throat Surgery |
| Indication | Adenoidectomy, Tonsillectomy |
| Surgical Step | Retraction |
| Type of Sterilization | not sterilized ex works |
| Data Matrix | Yes |
| Disposable / | |
| Reusable | reusable |
| Sterility | non-sterile |
| Product contains | |
| Latex | No |

По вопросам продаж и поддержки обращайтесь:

| | | | |
|-----------------------------|---------------------------------|--------------------------------|---------------------------|
| Алматы (7273)495-231 | Казань (843)206-01-48 | Новокузнецк (3843)20-46-81 | Смоленск (4812)29-41-54 |
| Архангельск (8182)63-90-72 | Калининград (4012)72-03-81 | Новосибирск (383)227-86-73 | Сочи (862)225-72-31 |
| Астрахань (8512)99-46-04 | Калуга (4842)92-23-67 | Омск (3812)21-46-40 | Ставрополь (8652)20-65-13 |
| Барнаул (3852)73-04-60 | Кемерово (3842)65-04-62 | Орел (4862)44-53-42 | Сургут (3462)77-98-35 |
| Белгород (4722)40-23-64 | Киров (8332)68-02-04 | Оренбург (3532)37-68-04 | Тверь (4822)63-31-35 |
| Брянск (4832)59-03-52 | Краснодар (861)203-40-90 | Пенза (8412)22-31-16 | Томск (3822)98-41-53 |
| Владивосток (423)249-28-31 | Красноярск (391)204-63-61 | Пермь (342)205-81-47 | Тула (4872)74-02-29 |
| Волгоград (844)278-03-48 | Курск (4712)77-13-04 | Ростов-на-Дону (863)308-18-15 | Тюмень (3452)66-21-18 |
| Вологда (8172)26-41-59 | Липецк (4742)52-20-81 | Рязань (4912)46-61-64 | Ульяновск (8422)24-23-59 |
| Воронеж (473)204-51-73 | Магнитогорск (3519)55-03-13 | Самара (846)206-03-16 | Уфа (347)229-48-12 |
| Екатеринбург (343)384-55-89 | Москва (495)268-04-70 | Санкт-Петербург (812)309-46-40 | Хабаровск (4212)92-98-04 |
| Иваново (4932)77-34-06 | Мурманск (8152)59-64-93 | Саратов (845)249-38-78 | Челябинск (351)202-03-61 |
| Ижевск (3412)26-03-58 | Набережные Челны (8552)20-53-41 | Севастополь (8692)22-31-93 | Череповец (8202)49-02-64 |
| Иркутск (395)279-98-46 | Нижний Новгород (831)429-08-12 | Симферополь (3652)67-13-56 | Ярославль (4852)69-52-93 |
| Россия (495)268-04-70 | Киргизия (996)312-96-26-47 | Казахстан (7172)727-132 | |