

# FRANKFURT

## Технические характеристики

### По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	



### BW080R

FRANKFURT Thyroidectomy Frame, detachable, non-sterile, reusable

Product search type	Thyroidectomy Frame
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
To be used with (code)	BW094R, BW093R, BW085R, BW090R, BW091R, BW092R, BW081R, BW086R, BW087R, BW088R
Connection	detachable
Discipline	Pediatric Surgery, General Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW081R

FRANKFURT Connecting Clamp, non-sterile, reusable

Product search type	Connecting Clamp
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
To be used with (code)	BW088R, BW087R, BW086R, BW093R, BW092R, BW091R, BW090R, BW081R, BW094R, BW085R, BW082R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW082R

FRANKFURT Frame For Children, frame, oval, depth: 180 mm, width: 240 mm, non-sterile, reusable

Product search type	Frame For Children
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	frame
Shape / Form	oval
Width (mm)	240.000 mm
Depth (mm)	180.000 mm
To be used with (code)	BV829R, BW081R, BW085R, BW093R, BV842R, BW090R, BW091R, BW092R, BW094R, BW086R, BW088R, BW087R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW085R

FRANKFURT Retractor, blade, 180 mm, depth: 14 mm, width: 17 mm, non-sterile, reusable



Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	180.00 mm
Width (mm)	17.000 mm
Depth (mm)	14.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW086R

FRANKFURT Retractor, blade, 180 mm, depth: 17 mm, width: 16 mm, non-sterile, reusable



Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	180.00 mm
Width (mm)	16.000 mm
Depth (mm)	17.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW087R

FRANKFURT Retractor, blade, 180 mm, depth: 18 mm, width: 22 mm, non-sterile, reusable




Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	180.00 mm
Width (mm)	22.000 mm
Depth (mm)	18.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW088R

FRANKFURT Retractor, blade, 180 mm, depth: 30 mm, width: 19 mm, non-sterile, reusable




Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	180.00 mm
Width (mm)	19.000 mm
Depth (mm)	30.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW090R

FRANKFURT Retractor, blade, 140 mm, depth: 30 mm, width: 20 mm, non-sterile, reusable




Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	140.00 mm
Width (mm)	20.000 mm
Depth (mm)	30.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW091R

FRANKFURT Retractor, blade, 120 mm, depth: 30 mm, width: 40 mm, non-sterile, reusable



Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	120.00 mm
Width (mm)	40.000 mm
Depth (mm)	30.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW092R

FRANKFURT Retractor, blade, 180 mm, depth: 35 mm, width: 11 mm, non-sterile, reusable



Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	180.00 mm
Width (mm)	11.000 mm
Depth (mm)	35.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	General Surgery, Pediatric Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW093R

FRANKFURT Retractor, blade, 140 mm, depth: 55 mm, width: 20 mm, non-sterile, reusable



Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	140.00 mm
Width (mm)	20.000 mm
Depth (mm)	55.000 mm
To be used with (code)	BW082R, BW081R, BW080R
Discipline	Pediatric Surgery, General Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### BW094R

FRANKFURT Retractor, blade, 120 mm, depth: 55 mm, width: 30 mm, non-sterile, reusable



Product search type	Retractor
Author	FRANKFURT
Aesculap Product	
Groups	General Surgical Instruments
Type	blade
Length	120.00 mm
Width (mm)	30.000 mm
Depth (mm)	55.000 mm
To be used with (code)	BW080R, BW081R, BW082R
Discipline	Pediatric Surgery, General Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



## BW099R

Fixation Element, OR-table fixation, non-sterile, reusable

---

Product search type	Fixation Element
Aesculap Product	
Groups	General Surgical Instruments
Type	OR-table fixation
Discipline	Pediatric Surgery, General Surgery
Surgical Step	Retraction
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No

**По вопросам продаж и поддержки обращайтесь:**

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	