

Steinmann

Технические характеристики

По вопросам продаж и поддержки обращайтесь:


Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	

Steinmann pins are available in diameter 4,0 , 4,5 , and 5,0 mm. All diameters are available in the length 100 up to 250 mm.



KX108S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 80 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	80.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX110S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 100 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	100.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX112S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 120 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	120.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX114S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 140 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	140.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX115S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 150 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	150.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX116S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 160 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	160.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX118S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 180 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	180.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable




KX120S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 200 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	200.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable




KX121S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 210 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	210.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX125S

STEINMANN Pin, trocar tip, triangular, 4 mm diam., 250 mm



Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	250.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX225S

STEINMANN PIN 4.5X150MM TROCAR PT.TRIANG

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	150.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular



KX228S

STEINMANN PIN 4.5X180MM TROCAR PT.TRIANG

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	180.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular



KX230S

STEINMANN PIN 4.5X200MM TROCAR PT.TRIANG

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	200.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular



KX235S


STEINMANN PIN 4.5X250MM TROCAR PT.TRIANG

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	250.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular



KX245S

STEINMANN Pin, trocar tip, triangular, 5 mm diam., 150 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	150.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX248S

STEINMANN Pin, trocar tip, triangular, 5 mm diam., 180 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	180.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX250S

STEINMANN Pin, trocar tip, triangular, 5 mm diam., 200 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	200.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable




KX255S

STEINMANN Pin, trocar tip, triangular, 5 mm diam., 250 mm




Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	250.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX260S

STEINMANN Pin, trocar tip, triangular, 5 mm diam., 300 mm



Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	General Surgical Instruments
Shape / Form	triangular
Length	300.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Surgical Step	Dissection-Bone
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable



KX261

Protection cap for pin, plastic, 4 mm diam.

Product search type	Pin
Aesculap Product	
Groups	General Surgical Instruments
Diameter (mm)	4.00 mm
Material	plastic
Surgical Step	Dissection-Bone
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	reusable



KX262

Protection cap for pin, plastic, 5 mm diam.

Product search type	Pin
Aesculap Product	
Groups	General Surgical Instruments
Diameter (mm)	5.00 mm
Material	plastic
Surgical Step	Dissection-Bone
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	reusable



KX308S

STEINMANN PIN 4.0X180MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	180.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX310S

STEINMANN PIN 4.0X200MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	200.00 mm
Diameter (mm)	4.00 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX325S

STEINMANN PIN 4.5X150MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	150.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX328S

STEINMANN PIN 4.5X180MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	180.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX330S

STEINMANN PIN 4.5X200MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	200.00 mm
Diameter (mm)	4.50 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX345S

STEINMANN PIN 5.0X150MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	150.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX348S

STEINMANN PIN 5.0X180MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	180.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX350S

STEINMANN PIN 5.0X200MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	200.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



KX360S

STEINMANN PIN 5.0X250MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	250.00 mm
Diameter (mm)	5.00 mm
Tip profile	trocar tip
Material	stainless implant steel
Data Matrix	No
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



LX757S

STEINMANN PIN 1.6X55MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	55.00 mm
Diameter (mm)	1.60 mm
Tip profile	trocar tip
Material	stainless implant steel
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread



LX758S

STEINMANN PIN 1.6X65MM W/THREAD TRIANG.

Product search type	Pin
Author	STEINMANN
Aesculap Product	
Groups	Osteosynthesis
Length	65.00 mm
Diameter (mm)	1.60 mm
Tip profile	trocar tip
Material	stainless implant steel
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	disposable
Tool hub	triangular
Type	with Centre Thread

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	