

# Технические характеристики

## По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	



### EJ126R


Trocar, 6 mm, 4,5 mm, 3,3 mm, 1,9 mm diam., set of 4 tips, one handle

Product search type	Trocar
Aesculap Product	
Groups	General Surgical Instruments
Type	with 4 tips
Diameter (mm)	6.00 mm, 4.50 mm, 1.90 mm, 3.30 mm
Tip profile	triangular tip
Configuration	Set of 4 tips, one handle
Surgical Step	Others
Discipline	Ear-Nose-Throat Surgery, General Surgery, Maxillo-Facial Surgery
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable /	
Reusable	reusable
Sterility	non-sterile
Product contains	
Latex	No



### EJ305R

SINUSCOPY TROCAR 5MM 85MM



Product search type	Trocar
Aesculap Product	
Groups	General Surgical Instruments
Working length (mm)	85.0 mm
Working length (Inch)	3 1/4
Diameter (mm)	5.00 mm
Tip profile	triangular tip
Discipline	Ear-Nose-Throat Surgery, Maxillo-Facial Surgery
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable /	
Reusable	reusable
Sterility	non-sterile
Product contains	
Latex	No



### EJ315R

#### TILLEY LICHTWITZ ANTRUM TROCAR D:5 FG



Product search type	Antrum Trocar
Author	TILLEY LICHTWITZ
Aesculap Product	
Groups	General Surgical Instruments, ENT
Length	185.00 mm
Length (Inch)	7 1/4
Working length (mm)	100.0 mm
Working length (Inch)	4
Diameter (mm)	1.70 mm
Surgical Step	Others
Indication	Sinus And Hypophysectomy
Discipline	Maxillo-Facial Surgery, Ear-Nose-Throat Surgery
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### EJ316R

#### TILLEY LICHTWITZ ANTRUM TROCAR D:6 FG



Product search type	Antrum Trocar
Author	TILLEY LICHTWITZ
Aesculap Product	
Groups	ENT, General Surgical Instruments
Length	185.00 mm
Length (Inch)	7 1/4
Working length (mm)	100.0 mm
Working length (Inch)	4
Diameter (mm)	2.00 mm
Surgical Step	Others
Indication	Sinus And Hypophysectomy
Discipline	Ear-Nose-Throat Surgery, Maxillo-Facial Surgery
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### EJ317R

#### TILLEY-LICHTWITZ ANTRUM TROCAR D:7 FG



Product search type	Antrum Trocar
Author	TILLEY LICHTWITZ
Aesculap Product	
Groups	ENT, General Surgical Instruments
Length	185.00 mm
Length (Inch)	7 1/4
Working length (mm)	100.0 mm
Working length (Inch)	4
Diameter (mm)	2.30 mm
Surgical Step	Others
Indication	Sinus And Hypophysectomy
Discipline	Ear-Nose-Throat Surgery, Maxillo-Facial Surgery
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### EJ318R

#### TILLEY-LICHTWITZ ANTRUM TROCAR D:8 FG



Product search type	Antrum Trocar
Author	TILLEY LICHTWITZ
Aesculap Product	
Groups	ENT, General Surgical Instruments
Length	185.00 mm
Length (Inch)	7 1/4
Working length (mm)	100.0 mm
Working length (Inch)	4
Diameter (mm)	2.70 mm
Surgical Step	Others
Indication	Sinus And Hypophysectomy
Discipline	Maxillo-Facial Surgery, Ear-Nose-Throat Surgery
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No



### EK074R

#### TROCAR PIN BLUNT (HASSON) 12/110MM

Product search type	Trocar Pin
Aesculap Product	
Groups	Endoscopic Technology
Product family	Access/Trocars
Altern.	
Productname	Obturator
Working length (mm)	110.0 mm
Working length (Inch)	4 1/4
Diameter (mm)	12.00 mm
Tip profile	blunt (Hasson)
threaded/smooth	smooth
Colour	yellow
Material	stainless steel
To be used with (code)	EK034R
Minimally Invasive	
Surg. (MIS)	Yes
Surgical Step	Minimal-invasive access
Indication	Laparoscopy
Type of Sterilization	not sterilized ex works
Sales unit	package of 1 piece
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains	
Latex	No



### EK099R

#### CONE FOR 12MM TROCARS

Product search type	Cone Attachment
Aesculap Product	
Groups	Endoscopic Technology
Product family	Access/Trocars
Diameter (mm)	12.00 mm
threaded/smooth	threaded
Material	stainless steel
To be used with (code)	EK038R, EK034R
Minimally Invasive	
Surg. (MIS)	Yes
Surgical Step	Minimal-invasive access
Indication	Laparoscopy
Type of Sterilization	not sterilized ex works
Data Matrix	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains	
Latex	No

**По вопросам продаж и поддержки обращайтесь:**

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	